1 OBJECTIVE

To determine the Young's modulus of material of given bar by Koenigs method.

2 REQUIREMENTS

Bar (metal/wooden), knife edges, mirrors, scale, telescope

3 INTRODUCTION




Figure 1: The experimental setup

The experimental arrangement is as shown in the figure. The bar is mounted on a knife edge which carries the load on a pan P. At the ends of the bar two mirrors M_1 and M_2 , almost normal to the bar, but slightly displaced, are fixed to enable a scale S to be seen in the telescope T, the light from S having suffered 2 reflections.

The telescope carries a crosswire in the eye-piece and the apparatus is arranged so that scale divisions as seen in the telescope coincides with crosswire. If now the bar is loaded with say 0.5 or 1 kg, as a result of the depression produced, the scale division viewed in the telescope will be altered by the particular divisions. One can determine the Young's modulus using the relation:

$$Y = \frac{3WL^2(2D + \alpha)}{2bd^3x} \tag{1}$$

Koenigs Method Semester 6

where,

W = mg where m is load in kg

 $l \rightarrow \text{distance between the knife edges in m}$

 $D \rightarrow \text{distance between scale}$ and the more remote mirror M_2 in m

 $\alpha \rightarrow$ distance between the mirrors in m

 $b \rightarrow \text{breadth of the bar in m}$

 $d \rightarrow$ Thickness of the bar in m

4 PROCEDURE

1. Arrange the apparatus as shown in figure.

- 2. Adjust the telescope to get a clear image of the scale on the mirror M_1 .
- 3. Place the 0.5 kg load in the pan bar and take the scale readings which coincides with crosswire. Note the value.
- 4. Repeat the above step by adding equal loads upto 3.5 kg and note the corresponding depression.
- 5. Find the values of l, D and α using a meter scale, and values of b and d using a screw gauge.
- 6. Plot a graph of load versus depression. The slope of the graph gives the value of $\left(\frac{m}{r}\right)$. Substituting it in the equation, Y can be obtained.

5 OBSERVATIONS

		Reading of scale in mirror (cm)			
Sl.no.	Load				Depression
	$\mid m \mid$	Load increase	Load decrease	Average	x
	(kg)				(cm)


Distance between the knife edged, $I = \underline{\hspace{1cm}}$
Distance between scale and the mirror, $M_2 =$
Distance between mirrors, $\alpha = \underline{\hspace{1cm}}$
Breadth of the bar, $b = \underline{\hspace{1cm}}$
Thickness of the bar, $d = \underline{\hspace{1cm}}$

Koenigs Method Semester 6

6 CALCULATIONS

The plot of load versus depression gives a straight line. The slope of this line is given by

slope =
$$\frac{\Delta x}{\Delta m}$$



From equation (1),

$$Y = \frac{3Wl^2(2D + \alpha)}{2bd^3x} = \frac{3mgl^2(2D + \alpha)}{2bd^3x}$$

$$\therefore Y = \frac{3gl^2(2D + \alpha)}{2bd^3x} \left(\frac{m}{x}\right)$$

$$\Rightarrow Y = \frac{3gl^2(2D + \alpha)}{2bd^3x} \left(\frac{1}{\text{slope}}\right)$$

7 RESULT

The Young's modulus of the bar is _______Pa